



"People you can Trust"

ALLIED
MEDICAL



ALLIED
MEDICAL

ALLIED MEDICAL LIMITED



DR. VINOD KOHLI
(MD and Founder)

1982



2 Doctors
1 Client
Unlimited Passion

Today



500+ Professionals
36000+ Domestic Clients
300+ International Clients
Unlimited Passion

About Us – “who we are, what we do”

Since its inception in 1982, Allied Meditec under the guidance of its founder Dr. Vinod Kohli (MD & Founder) an eminent Consultant - Anaesthesia & Intensive Care with an experience of over 50+ years in the field has marked its presence in the Healthcare Industry as technically most advanced manufacturer of Anaesthesia, Operation Theatre, Intensive Care, Emergency equipment and other lifesaving medical devices.



Allied Meditec was established with an aim to deliver advanced yet cost effective medical solutions. Allied Meditec is an ISO9001 and ISO13485 accredited company for last 20+ years and majority of the products manufactured are certified with European CE International quality standards. In addition to international certifications Allied Meditec also has Indian CDSCO registration and BIS approval for its products.

Global R&D Centres


At Allied Meditec, Research and Development (R&D) division drives our industry positioning as a leading medical manufacturer for Global Market pushing boundaries of science to discover and develop breakthrough solutions in the medical field. We help to shift the workload of user authentication from human being to technology as we are continuously creating an innovative set of technologies. Allied Meditec has three research laboratories located at London (United Kingdom), Rajasthan, Haryana, and upcoming centre at Uttar Pradesh in India. These centres work on new and existing product development making it easier for medical fraternity to have more intuitive approach for patient treatment.



Self Reliant Approach

Allied Meditec's core competency resides in its centralized ownership of every operational facet, spanning from component manufacturing through the conclusive stages of product development. The entire continuum, encompassing design and development, component manufacturing workshop, assembly, and quality control, is seamlessly executed in-house.

MILESTONES & AWARDS



The Most Promising Company of the Year 2022-23 in Medical Devices – Medgate Today



MSME of the Year 2022-23 (Manufacturing) Haryana



Udyog Award for Excellence in Manufacturing of Ventilators - Surgical Manufacturers & Traders Association (SMTA)



Gold Award in the Medium Enterprise Category CII-AMTZ MedTech Quality Champion Awards 2024



Lifetime Achievement Award for the Year 2025 – Medgate Today



Sales & Admin, Gurgaon



Assembly Line, Karoli



Training Centre,
Karoli



ANAESTHESIA



Allied Meditec Neptune Prime

- Anaesthesia Workstation

- ✦ 3 gas anaesthesia system
- ✦ Multiple ventilation modes - V-CMV, P-CMV, V-SIMV, P-SIMV, PCV-VG, PSV & Manual/Spont
- ✦ Electronic flowmeters with electronic gas mixing
- ✦ 15.6" TFT touch screen

Allied Meditec Neptune Plus

- Anaesthesia Workstation

- ✦ 3 gas anaesthesia system
- ✦ Multiple ventilation modes - VCV, PCV, SIMV-VC, SIMV-PC, CPAP/PS, PCV-VG, Man/ Spont
- ✦ Twin vaporiser mounting
- ✦ 12.1 inches capacitive touch screen
- ✦ Digital Flowmeter with electronic gas mixing (optional)

Jupiter Plus

- Anaesthesia Workstation

- ✦ 3 Gas anaesthesia systems
- ✦ Multiple ventilation modes - VCV, PCV, P-SIMV+PSV, V-SIMV+PSV, PSV and Manual
- ✦ Integrated ventilator with 7" TFT color touch screen display
- ✦ Medical grade ABS plastic frame
- ✦ Work surface lighting

Jupiter

- Anaesthesia Workstation

- ✦ 3 gas anaesthesia system
- ✦ Built in Anti-hypoxia device
- ✦ 7 inches color display with capacitive touch screen
- ✦ Suitable for Adult and Pediatric use
- ✦ Settable tidal volume
- ✦ Medical grade ABS plastic framework



Jupiter Series

- Anaesthesia Machine

- ✦ 3 different variant to choose from
- ✦ Anti-hypoxic devices in all variants
- ✦ Robust framework
- ✦ Time tested design



Allied Meditec Pisces

- Anaesthesia Vaporiser

- ✦ Among top 5 Anaesthesia Vaporisers in the World
- ✦ Available for Isoflurane and Sevoflurane
- ✦ Single hand operation
- ✦ Multiple mounting & filling options
- ✦ Interlock mechanism

Meditec ME-SORB™

- Carbon Dioxide Absorbent

- ✦ Special "D" profiles of the particles provides high ratio of surface area to volume.
- ✦ Low dust. Less than 0.3%
- ✦ High CO₂ absorption resulting low cost
- ✦ Low Interaction with the agent



Allied Meditec 2800

- Anaesthesia Ventilator


- ✦ Pneumatically driven ventilator
- ✦ Audio visual safety alarms
- ✦ Integrated rechargeable battery
- ✦ VCV, PCV, V-SIMV+PS, P-SIMV+PS, PSV
- ✦ High-luminance 7" TFT display screen



Allied Meditec Medlite

- Portable Anaesthesia Machine

- ✦ Ideal for OR, veterinary clinics, ambulance, military use, disaster management & emergencies
- ✦ Specially designed portable carry case
- ✦ Selectatec vaporiser manifold
- ✦ Multiple gas options



SharpScope 100

- Video Laryngoscope

- ✦ 3.5" non-reflective touch screen.
- ✦ Tilt and Rotatable Display
- ✦ 16 GB internal memory with expansion upto 256 GB
- ✦ Compatible with Macintosh disposable blades
- ✦ 4 Hours battery backup



Allied Meditec 1700

- Ventilator

- ✦ Electronically driven ICU ventilator
- ✦ Modes - PACV, PSIMV, VACV, VSIMV, AUTO (Emergency Quick Start), Spont+CPAP, PRVC, PRVC+SIMV, O₂ Stream
- ✦ Upgradable to SpO₂, EtCO₂ and ABIPAP+
- ✦ Optional - 12.1" TFT graphical display unit
- ✦ Trends and loops display



Allied Meditec Vista

- Neonatal to Adult Ventilator

- ✦ 15" large touch screen display
- ✦ Invasive & Non-Invasive Ventilation
- ✦ Suitable for adult to neonate patients
- ✦ Modes - PACV, PSIMV, VACV, VSIMV, Spont., Apnea, PRVC, Bi-level, AwPRV, AutoVent®, TCPL-SIMV, PRVC-SIMV, CPR, O₂ stream®
- ✦ SHFV, DHFV (optional)
- ✦ CO₂ & SpO₂ monitoring (optional)
- ✦ Medical air compressor (optional)



Allied Meditec IP200

- Volumetric Infusion Pump

- ✦ 3.5" TFT touch screen
- ✦ 8 operating modes
- ✦ Compatible with all IV set brands
- ✦ 11 adjustable, occlusion pressure settings
- ✦ 2000 event logs



Allied Meditec SP200

- Syringe Infusion Pump

- ✦ Syringe driven pump with automatic Syringe detection
- ✦ KVO (Keep vein open)
- ✦ Anti bolus system
- ✦ Audible and Visual alarms
- ✦ Programmable automatic & Manual bolus
- ✦ Compatible with all brands of syringes



Allied HFT

- High-Flow Therapy Device (HFNC)

- ✦ Flow Rate Range: 2 – 80L/min.
- ✦ 5" Color LCD touchscreen display
- ✦ Built-in turbine with electronic blender
- ✦ O₂% range: 21% – 100%



Allied Meditec IC20

- Infusion Workstation with Control Panel

- ✦ Large display screen
- ✦ Suitable for mounting up to 15 infusion pumps
- ✦ Wi-Fi and network interface
- ✦ Realtime data monitoring

MONITORING & ECG

Allied Meditec M751

- Pre-configured Patient Monitor

- ✦ 12.1 inches LED touch screen
- ✦ Upgradable to IBP, EtCO₂ and AG
- ✦ 480 hours of graphical/ tabular trends
- ✦ Advanced monitoring options designed for ICU, CCU & OR



Allied Meditec M707

- Upgradable Patient Monitor


- ✦ 5.5 inches LED display that can be used as plug transport module for M771 & M797
- ✦ 7-parameters monitoring capabilities with EtCO₂ and IBP as optional
- ✦ Provision of an additional rack for integration of advanced parameter modules



Allied Meditec M797

- Modular Patient Monitor

- ✦ 19 inches LED display
- ✦ Defibrillation interface and Cautery protection
- ✦ Advance parameter module upgradability: Respiratory mechanics, Regional Oximetry, PiCCO



Allied Meditec M700CMS

- Central Monitoring System

- ✦ Wired or wireless networking
- ✦ Display 17", 19", 21" and 24" TFT screen
- ✦ Shows maximum 32 patients information
- ✦ Intellectualized Bi-directional communication



Allied Meditec M300 Series

- Vital Signs Monitor

- ✦ 8" high resolution screen
- ✦ Portable, light weight and sturdy design
- ✦ Flexible parameter configuration for different clinical environment
- ✦ ECG, SpO₂, NIBP, Resp. & Temp. Monitoring



Allied Meditec E417

- Electro-Cardiograph

- ✦ Colored 8.4" Touch Screen
- ✦ 12 leads ECG simultaneous acquisition
- ✦ Alphanumeric keypad
- ✦ Patient data export in PDF format



Allied Meditec E427

- Electro-Cardiograph

- ✦ Colored 5.6" LCD display
- ✦ Automatic data interpretation
- ✦ In-built thermal printer
- ✦ 12 leads ECG simultaneous acquisition

EMERGENCY



Allied Meditec M1700 - Transport Ventilator

- ✦ Electronically driven ICU ventilator
- ✦ Modes - PACV, PSIMV, VACV, VSIMV, AUTO (Emergency Quick Start), Spont+CPAP, PRVC, PRVC+SIMV, O₂ Stream
- ✦ Upgradable to SpO₂, EtCO₂ and ABIPAP+
- ✦ Optional - 12.1" TFT graphical display unit
- ✦ Trends and loops display



Cardiasafe-I - Biphasic Defibrillator Monitor

- ✦ 7 inches display with intuitive 1-2-3 operation
- ✦ ECG, Resp, Thermal recorder, 360 Joules
- ✦ Paddle contact indicator
- ✦ Sliding built-in paediatric paddles



Cardiasafe-II - Biphasic Defibrillator Monitor

- ✦ AED, Pacer, Resp, Thermal Recorder, 360J
- ✦ Battery backup up to 6hrs or 200 shocks (Defib mode)
- ✦ Optional Features: NIBP, EtCO₂, SpO₂
- ✦ Interfacing: USB, RJ45, Multi-functional connector



Cardiagard-I

- Automated External Defibrillator (AED)

- ✦ Shock Control : Semi-automatic
- ✦ Cardiagard-I : LED lights, graphic panel and voice prompt
- ✦ Battery : new battery can discharge 200 ± 10 times at 150J
- ✦ At least 8 hours of ECG and event/ incident recording



Allied Meditec 1100

- Portable Ventilator

- ✦ Pneumatically driven and electronically controlled
- ✦ 7 inch LCD display
- ✦ VCV, SIMV, A/C, PCV, MANUAL, CPAP/PSV, S/T, S-mode, T-mode
- ✦ Invasive and Non-Invasive mode



Cardiagard-II

- Automated External Defibrillator (AED)

- ✦ Shock Control : Semi-automatic
- ✦ Cardiagard-II : 5"Color LCD screen, animation, voice and bilingual text prompts
- ✦ Battery : new battery can discharge 200 ± 10 times at 150J
- ✦ At least 8 hours of ECG and event/ incident recording



Allied Meditec 1200

- Portable Resuscitator

- ✦ Pneumatically driven and electronically controlled
- ✦ IPPV ventilation mode with electronically controlled alarms
- ✦ Pressure manometer and venturi suction facility
- ✦ 60% & 100% FiO₂ settings



Medical Fair (Thailand)



Medcall (India)



Medtech Expo



CME (Rewari)




Arab Health (Dubai)



Medical Fair 2024 (India)



Medical Fair 2025 (India)



Africa Health (South Africa)

India - Central Government & Programs



DGAFMS



Border Security Forces



Central Reserve Police Force



Indo Tibetan Border Police



AIIMS



NIMHANS



IOCL



NLC



NEW FRONTIERS, NEW DREAMS



सेल SAIL

International Clients



Ministry of Interior, Egypt



Ministry of Internal Affairs, Uganda



Department of Health, Philippines



Al-Saadi Private Hospital, IRAQ



Shahabuddin Medical College and Hospital, Bangladesh



Nong Phai Hospital, Thailand



RSUD Pameungpeuk Indonesia



Khairkwa Medical Complex, Afghanistan



Chhinnamasta Hospital, Nepal



Charis Veterinary Clinic, Bahrain



Ministry of Health, Kenya



Ministry of Health, Kuwait



Ministry of Health, Bhutan

India - State Governments



Govt. of NCT, Delhi



CGMSC (State Government Undertaking)



Karnataka State Medical Supplies Corporation Limited



Uttar Pradesh University of Medical Sciences



Sree Chitra Tirunal Institute for Medical Sciences & Technology



India - Corporate Hospitals



Partners in Health



the art of healthcare



Sir Ganga Ram Hospital



Manipal Hospital INSPIRED BY LIFE



Caring for you... for life





ALLIED MEDICAL LIMITED

Corporate Office : 76, Udyog Vihar, Phase IV, Gurugram - 122 015, Haryana

Phone: +91 124 4111444, Mobile: +91 9667462828

Toll Free: 1800 102 7879, Email: sales@alliedmed.co.in

Plant – I : F-126-127, G133-134, Industrial Area, Karoli, Bhiwadi - 301707, Rajasthan

Plant – II : 77, Udyog Vihar, Phase-IV, Gurugram - 122015, Haryana

Plant – III : Medical Device Park, YEIDA, Jewar - 203209, Uttar Pradesh

Meditec International England Limited

8, Pinner View, Harrow, Middx, HA1 4QA, United Kingdom

Phone: +44 (0) 203 6423808, Mobile: +44 7551584162

Email: sales@meditecengland.co.uk, Website: www.meditecengland.co.uk

