


Transducers

C5-2XQ	C5-2Q	C5-1Q	C6-2MQ
			
2 - 5 MHz Curved Array	2 - 5 MHz Curved Array	1 - 5 MHz Curved Array	2 - 6 MHz 3D/4D Mechanical Curved Array
L12-5Q	L12-5HQ	L12-5WQ	L17-7HQ
			
5 - 12 MHz Linear Array	5 - 12 MHz Linear Array	5 - 12 MHz Linear Array	7 - 17 MHz High Frequency Linear Array
L17-7SQ	MC8-4Q	P5-1Q	P7-3Q
			
7 - 17 MHz High Frequency Compact Linear Array	4 - 8 MHz Micro - Convex Array	1 - 5 MHz Phased Array	3 - 7 MHz Pediatric Phased Array
E10-3Q	E10-3HQ	E10-3BQ	ECL12-3Q
			
3 - 10 MHz Endocavity Curved Array	3 - 10 MHz Endocavity Curved Array	3 - 10 MHz Endocavity Curved Array	3 - 12 MHz Endocavity Bi-Plane Array

Being Productive Being Trust

Acclarix[™] LX3
Diagnostic Ultrasound System





Stunning Clarity

Compacted with innovative Edan TAI technology and multiple imaging processing technologies, Acclarix LX3 could perfectly display ultrasound image in different modes, assisting sonographer to make more precise diagnosis.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
 - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
 - In Doppler, TAI automatically adjusts for flow state, providing improved continuity, border detection and filling of color



Fetal Face (eLive)




Fetal Heart




MSK



Splenic Hilum



TVD



Anatomic M

New Generation of Capability and Flexibility

Integrated with compact appearance, cutting-edge technology and intelligent workflow, Acclarix LX3 is to provide more valuable innovation. Its excellent image performance plays an important role in the prevention and early detection in primary medical care, enriching the diagnostic application and enhancing the clinical value.

The Virtue of Value

The Acclarix LX3 delivers unmatched value and performance across a broad range of applications:

- Abdomen
- Cardiac Screening
- Small Parts
- OB/GYN
- Vascular
- MSK
- Nerve
- Breast

Distinctive Design



Considerate gel warmer design for better experience



Customizable touch screen user panel to simplify diagnosis



5 active transducer ports to maximize clinical application



Keyboard with dedicated swivel for easy folding