

Hospital Beds

HI-TECH - 09 (B)

ICU Bed Electric (Abs Panels & Abs Railing)

- Frame work made of rectangular M.S Tube.
- ABS moulded head & foot panels.
- Telescopic IV Rod with two location.
- Overall Size : 210 L x 90 W x 60-80 H cms.
- Backrest, Knee Rest Trendlenburg / Reverse Trendlenburg and Hi-Low Positions obtained by Electric Actuator operated by hand remote system.
- 4 Section top made of perforated MS Sheet.
- ABS tuck away safety side railing.
- Bed mounted on 125 mm dia castor (two with brakes)
- Finish : Epoxy Powder Coated.

Optional

- CPR Function on backrest section
- Central Braking Castors
- ACP Control
- Manual Override

* Optional-Bed also available with weighing scale



HI-TECH - 10 (B)

ICU Bed Electric (Abs Panels & Collapsible Railing)

- Frame work made of rectangular M.S. Tube.
- Telescopic IV Rod with four locations.
- Backrest, Knee Rest Trendlenburg / Reverse Trendlenburg and Hi-Lo Position by Electric Actuator operated by hand remote system.
- Overall Size : 210 L x 90 W x 60 H cms.
- Collapsible type safety side railing : SS horizontal tubes aluminium alloy / fixed with brackets ; Easily detachable.
- Perforated CRC Sheet top in four sections.
- Bed mounted on 125 mm dia castors-two with brakes.
- ABS moulded head & foot panels.
- Finish : Epoxy Powder Coated.

Optional

- CPR Function on backrest section
- Central Locking
- ACP Control
- Manual Override

